


10L WATER STILL

Model: 406



FEATURES

- Water distillers are used to produce high quality water through the process of distillation
- This is an effective method of producing high quality water from a raw source and eliminating contaminants, minerals and trace elements
- For indoor use only
- Ambient temperature range +5°C to +40°C
- Relative humidity conditions: 80% RH for temperatures up to 31°C decreasing linearly to 50% RH at 40°C
- Safety cut out in the event of low water
- Easy access for cleaning and de-scaling tank and elements



NOTE

- 10ℓ/h water yield



ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO
YEAR WARRANTY ON COMPONENTS AND
DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.


www.scientific.co.za

10L WATER STILL

Model: 406

SPECIFICATIONS

| | |
|-----------------------|----------------------|
| POWER | 7500W |
| VOLTAGE | 380V / 50Hz |
| PACKAGING (W x D x H) | 385 x 385 x 925 (mm) |
| UNIT (W x D x H) | 350 x 270 x 830 (mm) |
| NET WEIGHT | 11kg |
| GROSS WEIGHT | 12kg |
| SHIPPING VOLUME | 0.14m ³ |
| CAPACITY | 10ℓ |
| WATER YIELD | 10ℓ/h |



NOTE

- When viewing the appliance from the front in its usual operating position, the width (W) of the product is the total distance from left to right; the depth (D) of the product is the total distance from the front to the back; the height (H) is the total distance from the bottom of the product to the top. Pictures may vary from actual product.

DISCLAIMER

- Output performance figures quoted are dependant on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference.

WARRANTY

- All Scientific equipment comes with a standard ONE year warranty. The 1 Year Extended Warranty is only valid for customers who have registered their warranty online at www.scientific.co.za within the first year after making their purchase.