

800ml MIXER

Model: 723



FEATURES

- The Scientific mixer is durable and robust for use in agricultural, microbiological and cement laboratories
- Rugged construction
- Small footprint

NOTE

- Supplied with a 800ml stainless steel mixing cup



ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO
YEAR WARRANTY ON COMPONENTS AND
DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.

www.scientific.co.za






800ml MIXER

Model: 723

SPECIFICATIONS

POWER	230V / 50Hz
VOLTAGE	350W
PACKAGING (W x D x H)	190 x 180 x 570 (mm)
UNIT (W x D x H)	180 x 260 x 470 (mm)
NET WEIGHT	4.5kg
GROSS WEIGHT	5.5kg
SHIPPING VOLUME	0.02m ³
CUP CAPACITY	800ml



NOTE

- When viewing the appliance from the front in its usual operating position, the width (W) of the product is the total distance from left to right; the depth (D) of the product is the total distance from the front to the back; the height (H) is the total distance from the bottom of the product to the top. Pictures may vary from actual product.

DISCLAIMER

- Output performance figures quoted are dependent on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference.

WARRANTY

- All Scientific equipment comes with a standard ONE year warranty. The 1 Year Extended Warranty is only valid for customers who have registered their warranty online at www.scientific.co.za within the first year after making their purchase.



**ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO
YEAR WARRANTY ON COMPONENTS AND
DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.**

www.scientific.co.za

Published Date: 2020.09.25