

## CERAMIC MAGNETIC STIRRER

Model: 725



## FEATURES

- ❑ Magnetic Stirrers are suitable for most types of laboratory applications where controlled stirring and heating is required
- ❑ Corrosion resistant ceramic heating plate
- ❑ Reliable heating elements
- ❑ Temperature range 50°C - 350°C
- ❑ Reliable solid state electronic speed control
- ❑ Stirring bar 40 x 8mm
- ❑ RPM: 100 - 2000
- ❑ Control of temperature via hotplate control or PT 100 inserted into the beaker

### NOTE

- ❑ Stand, clamp, external PT100 probe and ceramic stirrer bar included






ALL SCIENTECH EQUIPMENT COMES WITH A ONE  
YEAR WARRANTY ON COMPONENTS AND  
DEFECTIVE WORKMANSHIP

# CERAMIC MAGNETIC STIRRER

Model: 725

## SPECIFICATIONS

POWER	600W
VOLTAGE	230V / 50Hz
PACKAGING (W x D x H)	240 x 400 x 220 (mm)
UNIT (W x D x H)	190 x 350 x 400 (mm)
NET WEIGHT	8kg
GROSS WEIGHT	8.5kg
SHIPPING VOLUME	0.02m <sup>3</sup>
PLATE SIZE	190 x 190 (mm)
RPM	100 - 2000
LOAD CAPACITY	5kg



## TERMS & CONDITIONS

**NOTE:** Pictures may vary from actual product. | **DISCLAIMER:** Output performance figures quoted are dependent on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference. | **WARRANTY:** All Scientech equipment comes with a standard ONE year warranty.

REMEMBER TO COMPLETE THE ONLINE WARRANTY REGISTRATION



ALL SCIENTECH EQUIPMENT COMES WITH A ONE  
YEAR WARRANTY ON COMPONENTS AND  
DEFECTIVE WORKMANSHIP

Published Date: 2020.09.25